



ZIP-CoVx-P2 Test (SARS-CoV-2)

Key features:

- Molecular point-of-care test
- Direct detection of SARS-CoV-2 virus
- Two gene targets for test result confirmation
- Combined throat and nasal swab sampling
- Results within 7-30 minutes
- Sensitivity comparable to PCR
- Transport and storage at room temperature



Use with the ZIP-P2 Instrument

- Deployable at point-of-care settings
- Decentralised diagnosis
- Integrated sample preparation
- Screen-guided easy-to-use workflow with automated results and reporting



Product Specifications

Intended Use

The ZiP-CoVx-P2 test detects SARS-CoV-2 virus nucleic acids in combined throat and nasal swab samples to aid diagnosis of COVID-19 disease in symptomatic individuals, or to screen for SARS-CoV-2 infection in asymptomatic individuals. The test is intended for use in dedicated test spaces by healthcare professionals, or in pathology settings by laboratory trained professionals. The ZiP-CoVx-P2 test is performed using the ZiP-P2 instrument.

Sample Types

A combined oropharyngeal (throat) and bilateral mid-turbinate (nasal) sample, collected on a synthetic flocked swab. Dry swab samples must be used because swabs in liquid transport media may interfere with test performance.


Targets and Strains Detected

SARS-CoV-2; dual gene targets for *M* and *ORF1ab* genes. Detects at minimum, all variants of concern including: Alpha (B.1.1.7), Beta (B.1.351), Delta (B.1.617.2), and Omicron (BA.2, BA.4, BA.2.75, and XBB).

Endogenous control; human *RNaseP* gene. This is used to confirm sample adequacy and to ensure no test inhibition.

Test Workflow

Simple 4 step process, guided by on-screen instructions. Minimal training required.



Analytical and Clinical Performance

The overall analytical limit of detection (LoD) of the ZiP-CoVx-P2 test is 2,000 copies/swab (*M* gene 2,000 copies/swab, *ORF1ab* gene 4,000 copies/swab).

The overall clinical performance of the ZiP-CoVx-P2 test is 98.3% sensitivity (114/116) and 99.0% specificity (589/595) when compared to lab-based RT-PCR testing.

Interpretation of Results

Results are interpreted automatically by the ZiP-P2 instrument and shown on-screen.

Results are reported as POSITIVE++ (both SARS-CoV-2 gene targets are detected), POSITIVE + (one SARS-CoV-2 gene target is detected), NEGATIVE or INVALID.

This is a qualitative test and does not provide the quantitative value of the detected organism.

Storage and Stability

The ZiP-CoVx-P2 test is stable until its expiry date if stored between 2°C to 25°C.

Avoid direct light. Do not freeze.



THIS PRODUCT IS NOT AVAILABLE FOR PURCHASE BY THE GENERAL PUBLIC
THIS PRODUCT IS NOT COMMERCIALY AVAILABLE IN ALL MARKETS
ALWAYS FOLLOW THE DIRECTIONS FOR USE



ZIP Diagnostics Pty Ltd

24 Cromwell St, Collingwood, Victoria, 3066, Australia

Phone +61 (0) 3 8414 5770 | Email info@zipdiag.com

© 2023 ZiP Diagnostics Pty Ltd. All rights reserved.
D008366 v1.2 (Feb 2024)

