




ZIP-RPA Isothermal Amplification Kit

Strand-displacing and polymerisation enzymes and reagents engineered for rapid isothermal DNA amplification. Key features:

- Optimised for 37-42°C
- Results typically within 5-10 minutes
- Suitable for diagnostic settings requiring rapid results and low heating, or where elevated temperatures may be undesirable
- Use in next generation sequencing workflows

Ordering Information

- Standard supply is FOR RESEARCH USE ONLY
(Not for supply or use in US regions)
- Can migrate to a quality contract for supply and licensing for use in regulated animal health, food safety and human health applications
- Contact ZIP for order and supply scheduling
sales@zipdiag.com



ZIP Diagnostics Pty Ltd

24 Cromwell St, Collingwood, Victoria, 3066, Australia
Phone +61 (0) 3 8414 5770 | Email info@zipdiag.com

© 2023 ZIP Diagnostics Pty Ltd. All rights reserved.
D008560 v1.1 (Feb 2024)

